

**1. Part Number:** 5005A-250

**2. Photo:**



### 3. Product Information:

Part No.	5005A-250
Description	Stainless Steel Dipper – 250ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	250ml
Brimful Cup Volume:	250ml
Overall Length:	600mm
Max. Width (cup and handle)	73mm
Cup - Outer Diameter:	63mm
Cup - Height:	80mm
Nominal Weight:	555g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

### 4. Approval:

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 10-Aug-23