

1. Part Number: 5005B-100

2. Photo:



3. Product Information:

Part No.	5005B-100
Description	Stainless Steel Dipper – 100ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	100ml
Brimful Cup Volume:	100ml
Overall Length:	1000mm
Max. Width (cup and handle)	49mm
Cup - Outer Diameter:	39mm
Cup - Height:	100mm
Nominal Weight:	735g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 10-Aug-23