

1. Part Number: 5005B-50

2. Photo:



3. Product Information:

Description	Stainless Steel Dipper – 50ml – 1000mm
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	50ml
Brimful Cup Volume:	50ml
Overall Length:	1000mm
Max. Width (cup and handle)	42mm
Cup - Outer Diameter:	32mm
Cup - Height:	80mm
Nominal Weight:	710g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 10-Aug-23