


1. Part Number: 5005C-250

2. Photo:



3. Product Information:

Part No.	5005C-250
Description	Stainless Steel Dipper – 250ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	250ml
Brimful Cup Volume:	250ml
Overall Length:	1500mm
Max. Width (cup and handle)	73mm
Cup - Outer Diameter:	63mm
Cup - Height:	89mm
Nominal Weight:	1180g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 10-Aug-23