

**1. Part Number:** 8045P-1000S

**2. Photograph:**



**3. Product Information:**

**Description:**

SteriWare PharmaScoop® PP 1000ml - Sterile

**Use:**

Single use, disposable scoop

**Material of Construction:**

Polypropylene – rigid, opaque

- FDA acceptable CFR 177.1520
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant
- BSE/TSE free

**Method of Construction:**

Single piece, injection moulded

**Nominal Scoop Volume:**

1000ml

**Overall Length:**

324mm

**Handle Length:**

149mm

**Maximum Width:**

150mm

**Nominal Weight:**

138g

**Moulding Environment:**

Class 100,000 Cleanroom

**Primary Packaging Environment:**

Class 10,000 Cleanroom

**Individually Bagged?**

Yes (heat sealed PE/PET bag)

**Method of Sterilisation:**

Gamma Irradiation (25-45kGy)

**Number of Scoops per Box:**

25

**Recommended Storage Conditions:**

Dry and ambient temperature

**Shelf Life:**

5 years from the date of manufacture

**ATEX:**

This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 21-Nov-22