

1. Part Numbers: A826-12

2. Photograph:



3. Product Information:

Description:	Bucket 12L
Material of Construction:	316L stainless steel
Interior Finish:	Better than 1.0 microns Ra
Max. Diameter:	310mm
Height:	290mm
Graduations:	Every 1L (stamped)
Nominal Weight:	1.7kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-May-21