


**1. Part Number:** A866-1000

**2. Photograph:**



### 3. Product Information:

Description:	Beaker with Spout 1000ml
Method of Construction:	Crevise free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	1000ml
Brimful Volume:	1230ml
Dimensions:	
A - Internal Diameter:	98mm
B - Nominal Overall Diameter:	119mm
C - Nominal Height:	161mm
Graduation Marks:	Every 100ml from 400ml up to 1000ml
Nominal Weight:	520g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
 Position: Quality Manager  
 Date: 19-Oct-23