


1. Part Number: A866-2000

2. Photograph:




3. Product Information:

Description:	Beaker with Spout 2000ml
Method of Construction:	Crevise free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	2000ml
Brimful Volume:	2550ml
Dimensions:	
A - Internal Diameter:	140mm
B - Nominal Overall Diameter:	153mm
C - Nominal Height:	171mm
Graduation Marks:	Every 500ml from 1000ml up to 2000ml
Nominal Weight:	800g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 19-Oct-23