


1. Part Number: A866-200

2. Photograph:



3. Product Information:

Description:	Beaker with Spout 200ml
Method of Construction:	Crevice free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	200ml
Brimful Volume:	220ml
Dimensions:	
A - Internal Diameter:	59mm
B - Nominal Overall Diameter:	69mm
C - Nominal Height:	80mm
Graduation Marks:	None
Nominal Weight:	93g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 19-Oct-23