


1. Part Number: A866-3000

2. Photograph:



3. Product Information:

Description:	Beaker with Spout 3000ml
Method of Construction:	Crevice free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	3000ml
Brimful Volume:	3300ml
Dimensions:	
A - Internal Diameter:	140mm
B - Nominal Overall Diameter:	155mm
C - Nominal Height:	221mm
Graduation Marks:	Every 500ml from 1000ml up to 3000ml
Nominal Weight:	1000g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 19-Oct-23