


1. Part Number: A866-500

2. Photograph:




3. Product Information:

Description:	Beaker with Spout 500ml
Method of Construction:	Crevise free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	500ml
Brimful Volume:	610ml
Dimensions:	
A - Internal Diameter:	79mm
B - Nominal Overall Diameter:	99mm
C - Nominal Height:	121mm
Graduation Marks:	Every 100ml from 200ml up to 500ml
Nominal Weight:	315g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 19-Oct-23