

| | | | |
|---|--------------------------------|----------------------|------------|
|  | <h1>Product Specification</h1> | Specification Number | 8890P-500S |
| | | Revision | 01 |


Part No. 8890P-500S
Description: SteriWare® Powder Bottle 500ml
Material of Construction: Polypropylene (PP) – opaque

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Nominal Volume: 500ml
Brimful Volume: 540ml
Method of Construction: Single piece, injection moulded
Height: 106mm
Lid Diameter: 94mm
Nominal Weight: 83g
Assembly & Packing Environment: Class 100,000 Cleanroom
Individually Bagged? Yes (heat sealed PE/PET bag)
Method of Sterilisation: Gamma Irradiation (25-45kGy)
Number of containers per Box: 37
BSE/TSE Free: Yes
Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade mark belonging to Sampling Systems Ltd.

| | | |
|----------------------------|--|--|
| Specification Approved By: | Signature:  | Name: Jonathan Fenner Position: Director Date: 04 April 2017 |
|----------------------------|--|--|