



	<h1>Product Specification</h1>	Specification Number	A836-250
		Revision	01

Part No.	A836-250		
Description	Stainless Bottle 250 mL		
Nominal Volume:	250 mL		
Brimful Volume	380 mL		
Material of Construction:	Body:	316 stainless steel	
	Lid:	316 stainless steel	
	Clamp:	304 stainless steel	
	Gasket:	Silicone (conforms to FDA CFR 177.2600)	



Method of Construction:	Crevice free interior, welds ground & polished
Surface Finish:	Better than 0.5 microns Ra
Dimensions:	A = 64 mm B = 48 mm C = 97 mm D = 83 mm
Nominal Weight (body, lid, clamp & gasket):	710 g

Note: The neck of the container is fitted with a 2" diameter ferrule which conforms to BS4825-3

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 14-Sep-16
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	