



Metal Detectable

Reusable and Single use items for Laboratories, Quality Control and Production



World

Class Equipment

worldwide



Welcome to the New Metal Detectable Range catalogue.

After many requests from customers we are excited to introduce the new Metal Detectable range. It is an extensive range of x-ray and metal detectable products.

Every item in our metal detectable range is manufactured and packed in a Cleanroom to ensure product purity & quality.

Metal Detectable products are ideal for use in food, chemical and pharmaceutical industries.

Documentation can be easily downloaded from our website (www.sampling.com):

- Specification sheets
- Certificates of Conformity

If you can not see what you require, please contact us. We look at every application and do our very best to help you.

We look forward to hearing from you.

Jonathan Fenner

Simon Fenner



Made in the UK

Quality Guaranteed



The Sampling Systems Difference

Cleanroom Manufacture



100% Inspection



Individually Bagged



Full Documentation & Traceability



Over one million stock items



Fast delivery worldwide



Ultra Strong Packaging with Easy Opening



Highly Visible Colour



Specials and New Designs



Metal Detectable Scoop

High Quality Scoops for Dispensing and Sampling



Our highly visible blue Metal Detectable Scoops are ideal for food and pharmaceutical companies. The lightweight ergonomic design means they are easy to use. The distinctive blue colour ensures they can be easily seen should they be left in the powder (in line with HACCP guidelines).

- Moulded & packed in a Cleanroom for superior cleanliness
- Each scoop is individually wrapped in our ultra strong packaging that is designed to be easily opened
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming Polypropylene (PP)
- BSE/TSE free
- Large stocks for fast delivery worldwide
- DURABLE - can be used straight from the bag or repeatedly used



Part No.	Nominal Capacity	Colour	Material	Quantity per Pack
M1000B-30	30ml	Blue	PP	4
M1000B-60	60ml	Blue	PP	4
M1000B-125	125ml	Blue	PP	4
M1000B-250	250ml	Blue	PP	4
M1000B-500	500ml	Blue	PP	4
M1000B-1000	1000ml	Blue	PP	4
M1000B-2500	2500ml	Blue	PP	4

Metal Detectable Long Handle Scoop

Keeps hands out of the product

Long Handle Scoops are the ideal tool for reaching into deep containers. The extra long handle will help operators keep their hands out of the product. The other advantage of a long handle is that it makes the scoop more visible if it is left in the powder.

Metal Detectable Long Handle Scoops are manufactured and packed in a Cleanroom this ensures that the item is clean and ready for use.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 food grade conforming PP
- BSE/TSE Free
- Full Batch Traceability
- DURABLE - can be used straight from the bag or repeatedly used




New




Part No.	Nominal Scoop Volume	Overall Length	Overall Width	Material	Quantity per Pack
M1040B-100	100ml	400mm	50mm	PP	4


Download our latest **Brochures** at www.sampling.com



New Range
Brochure



New SteriWare
Brochure



New Stainless Ware
Brochure

Metal Detectable Spoon

Highly Visible Sampling Spoon



The Metal Detectable Spoon from Sampling Systems is ideal for sampling small amounts of powder and for dispensing small quantities.

The bright blue colour ensures that the spoons are highly visible and easy to spot if left in most powders or granules (in line with HACCP guidelines).

Each spoon is individually bagged and can be used straight from the bag.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 food grade conforming Polypropylene (PP)
- BSE/TSE Free
- Full Batch Traceability
- Available from stock



Part No.	Nominal Spoon Volume	Colour	Overall Length	Material	Quantity per Pack
M1006B-10	10ml	Blue	185mm	PP	4

Metal Detectable Micro-Spatula

Ideal for handling small quantities of powders & granules



The Metal Detectable Micro-Spatula is ideal for sampling small volumes. Its narrow sampling blade enables samples to be taken from very small containers such as vials.

Alternatively it can be used for dispensing small amounts of powders and granules.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Nominal Sample Volume	Colour	Material	Blade Width x Length	Overall Length	Quantity per Pack
M1013B-170	2ml	Blue	PP	13 x 40mm	170mm	4

Metal Detectable Volumetric Spoon

Accurate measuring spoons



Metal Detectable Volumetric Spoons allow small volumes to be accurately measured. The Volumetric Spoons are designed to be used straight from the bag. Suitable for use with powders and liquids.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU food grade 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability

NB. Volumetric spoons are not certified measuring devices. User validation is required. If in doubt please contact us



Part No.	Flat Fill Volume	Bowl Diameter	Overall Length	Material	Quantity per Pack
M1030B-25	2.5ml	23mm	160mm	PP	4

Metal Detectable Long Handled Spoon


Long Reach Spoon



Our Metal Detectable Long Handle Spoons are ideal when you need to reach into deep containers. The long handle and the bright colour help them to be quickly and easily seen if they are left behind in the powder.

The ergonomic handle is comfortable to use with the spoon slightly angled.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU food grade 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Flat Fill Volume	Handle Length	Material	Quantity per Pack
M1036B-20	20ml	300mm	PP	4

Metal Detectable Jug

Graduated for easy use

Our range of high quality Metal Detectable Jugs are designed to be used straight from the bag, no pre-cleaning is required. The specially designed ergonomic handle is rounded to make it comfortable to hold and use.

These jugs have a sturdy base and a pouring spout.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability
- Graduations are on the outside of the jug



Part No.	Nominal Volume	Diameter	Height	Material	Quantity per Pack
M1044B-500	500ml	90mm	130mm	PP	4
M1044B-1000	1000ml	120mm	170mm	PP	4

Metal Detectable Knife

Ideal for sample preparation

The Metal Detectable Knife can be used for preparing samples as well as general laboratory use.

Each knife is moulded and packed in a Cleanroom and can be used straight from the bag.

- Moulded & packed in a Cleanroom for superior cleanliness
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Individually Bagged	Overall Length	Blade Length	Material	Quantity per Pack
M1052B-100	Yes	160mm	60mm	PP	4

Metal Detectable Pallet Knife & Stirrer

Ideal for Production and Laboratory Use



The Metal Detectable Pallet Knife is very effective when working with creams and ointments, it is also a very efficient stirrer. The ergonomic handle makes them comfortable to use.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Blade Length	Blade Width	Overall Length	Material	Quantity per Pack
M1016B-200	200mm	45mm	300mm	PP	4

Metal Detectable Scraper

Cleanroom manufactured



The perfect tool for removing material from flat surfaces including tables, work tops and general equipment.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Blade Length	Blade Width	Overall Length	Material	Quantity per Pack
M1019B-110	110mm	110mm	190mm	Blue PP	4

DISCOUNTS AVAILABLE FOR BULK ORDERS - PLEASE CONTACT US

Metal Detectable Powder Funnel

Wide spout for efficient powder flow

The Metal Detectable Powder Funnel is ideal for transferring powders and granules. The wide spout helps to prevent powder blockages. Metal Detectable Powder Funnels are made from chemically resistant food grade PP.

- Individually bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- Moulded & packed in a Cleanroom for superior cleanliness
- Available from stock for fast delivery worldwide



New



Part No.	Top Diameter	Spout Diameter	Material	Quantity per Pack
M1051B-200	200mm	40mm	PP	4

Metal Detectable Liquid Funnel

Safely transfer liquids and suspensions

Our Metal Detectable Liquid Funnels are ideal when decanting liquids and suspensions. The Funnels have a small side handle to keep hands away from the product contact area when the funnels are in use.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



New



Part No.	Top Diameter	Spout Diameter	Overall Height	Material	Quantity per Pack
M1023B-100	100mm	12mm	125mm	PP	4

Metal Detectable Ladle

Efficient surface sampling of liquids

The Metal Detectable Ladle is designed to quickly and efficiently take surface samples from liquids. They are moulded and packed in a Cleanroom and so can be used straight from the bag.

The one piece moulded design has been specially made to fit through a standard 2" lid bung. 2" bungs are standard features on many barrels and kegs.



- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE free
- Full Batch Traceability
- Available from stock



Part No.	Nominal Volume	Diameter	Colour	Material of Construction	Quantity per Pack
M1033B-50	50ml	48mm	Blue	PP	4

Metal Detectable Bottle

Individually Bagged Bottles Made in a Cleanroom

Metal Detectable Bottles are designed to meet the requirements of the pharmaceutical, food and other hygienic industries.

The bottles are moulded and packed in a Class 100,000 Cleanroom. They have a wide opening and so are ideal when working with powders.

- Powder and splash tight
- Moulded & packed in a Cleanroom for superior cleanliness
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability
- Individually bagged



Part No.	Nominal Volume	Overall Height	Outer Diameter	Material	Quantity per Pack
M1047B-250	250ml	88mm	76mm	PP	4
M1047B-500	500ml	106mm	94mm	PP	4

Metal Detectable V-Spatula


New

For Perfect Sampling



Quick surface and core sampling is easy with our V-Spatula. Available in 2 lengths to offer shallow or deep sampling. They can also be used to take samples from sacks.

- Moulded & packed in a Cleanroom for superior cleanliness
- Ideal for general sampling and dispensing
- Suitable for free flowing and cohesive powders
- Ergonomic handle
- Free standing design



Part No.	Blade Length	Blade Width	Overall Length	Colour	Material	Quantity per Pack
M1010B-150	150mm	20mm	230mm	Blue	PP	4
M1010B-250	250mm	20mm	330mm	Blue	PP	4

Address:
Sampling Systems Ltd.
2 Trillennium Park
Gorsey Lane
Coleshill
B46 1JU
UK

Web: sampling.com
Tel: +44 (0)1675 466992
Email: info@sampling.com

Notes:
All goods are sold subject to our normal Terms & Conditions. Contact us for a copy.

