


1. Part Number: 1030BT all sizes

2. Product Information:



Part No:	Nominal Volume	Weight
1030BT-2	2ml	190g
1030BT-3	3ml	180g
1030BT-5	5ml	168g
1030BT-8	8ml	115g
1030BT-10	10ml	128g

Description

Tips for use with 19mm diameter PowderThief Bodies

Materials of Construction:

316L stainless steel

Finish:

<1 microns Ra

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin

Recommended Storage Conditions:

Dry and ambient temperature

3. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Aug-23