


**1. Part Number:** 1030CT all sizes

**2. Product Information:**



Part No:	Nominal Volume	Weight
1030CT-10	10ml	299g
1030CT-20	20ml	221g
1030CT-25	25ml	182g

**Description:** Tips for use with 25mm diameter PowderThief Bodies

**Materials of Construction:** 316L stainless steel

**Finish:** <1 microns Ra

**BSE/TSE statement:** All polishing compounds used in the manufacture of this product are of vegetable origin

**Recommended Storage Conditions:** Dry and ambient temperature

**3. Approval:**

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 29-Jan-21