

**1. Part Number:** 1035E-1500

**2. Photo:**



### 3. Product Information:

Description	Sample Sleeve Sampler Handle Extension Rod 1500mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Overall Length:	1525mm
Weight:	950g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature
Note:	This Extension Rod is designed to be used with either the 1000mm handle (1035A-1000) or 1500mm handle (1035A-1500)

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Nov-24