

1. Part Number: 1074A-350

2. Photograph:



3. Product Information:

Description	MicroThief 350mm
Materials of Construction:	- 316L stainless steel - Nylon Clutch (non product contact part)
Nominal Sample Volume:	Variable between 0.01 and 0.1ml
Overall Length:	450mm
Length of Sampling Tube:	350mm
Diameter of Sampling Tube:	5mm
Weight:	185g
Finish:	<0.8 microns Ra
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24