


1. Part Number: 1210B-1000

2. Photo:



3. Product Information:

| | |
|-------------------------------|----------------------------|
| Description | Slot Sampler 1000mm |
| Materials: | 316L stainless steel |
| External Polish: | Better than 0.8 microns Ra |
| Internal Polish: | None |
| Outer Diameter of Outer Tube: | 25mm |
| Inner Diameter of Inner Tube: | 17mm |
| Overall Length: | 1025mm |
| Slot Length: | 120mm |
| Slot Width: | 15mm |
| Number of Slots | 4 |
| Slot Positions: | |



| | |
|---------------------------------|---|
| Nominal Weight of Sampler: | 1.9kg |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this sampler are free of animal products |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 10-May-23