


**1. Part Number:** 1210B-1500

**2. Photo:**



### 3. Product Information:

Description	Slot Sampler 1500mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Internal Polish:	None
Outer Diameter of Outer Tube:	25mm
Inner Diameter of Inner Tube:	17mm
Overall Length:	1525mm
Slot Length:	120mm
Slot Width:	15mm
Number of Slots	5
Slot Positions:	



Nominal Weight of Sampler:	2.9kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 20-Mar-20