


1. Part Number: 1220R-1

2. Photo:



3. Product Information:

Description

Emptying Rack for Pocket Sampler 800mm

Material of Construction:

Funnels:	316L stainless steel
Bar:	316L stainless steel
Stands:	316L stainless steel
End Stop:	316L stainless steel
Thumb Screw:	Stainless Steel
Bottle 250ml:	Polypropylene

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004
- BSE/TSE free

Finish for Product Contact Parts:

Better than 0.5 microns Ra

Overall Length:

800mm

Max. Height:

240mm

Max. Width:

180mm

Max. Height Under Funnel:

105mm

Number of Funnels

3

Number of Bottles:

3

Nominal Weight of Rack:

3.1 Kg

Recommended Storage Conditions:

Dry and ambient temperature

4. Spare Parts:

Bottle 250ml

Part Number: 8185P-250 (Supplied as Box of 100)

5. Compatibility:

Pocket Samplers less than 1m in length = Buy 1 rack

Pocket Samplers greater than 1m in length = Buy 2 racks

A single Rack is suitable for all Pocket Samplers up to 1m in length (except 1220B-300). Two racks will be required for use with Pocket Samplers longer than 1m.

NOT suitable for use with: 1220B-300

Pocket Sampler Part Number	Sampler Length	Sampler Diameter	Number of Emptying Racks Required
1220A-600	600mm	25mm	1
1220A-900	900mm	25mm	1
1220A-1500	1500mm	25mm	2
1220B-300	300mm	25mm	Not Suitable
1220B-600	600mm	25mm	1
1220B-900	900mm	25mm	1
1220B-1500	1500mm	25mm	2
1220J-900	900mm	40mm	1
1220J-1500	1500mm	40mm	2
1220J-2000	2000mm	40mm	2

6. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24