


1. Part Number: 1360A-1003

2. Photo:



3. Product Information:

Description: AccuPod® Sample Thief 1000mm, 3 ports
Material of Construction: 316L stainless steel
Method of Construction: Welds ground and polished
Surface Finish: Better than 0.8 microns Ra (tube exteriors)
Overall Length: 1000mm
Outer Tube Diameter: 25mm
Number of Sampling Ports: 3
Position of Sampling Ports:



Nominal Weight: 2.5 kg
BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions: Dry and ambient temperature
First Use: This sampler should be thoroughly cleaned before being used for the first time
Note: This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24