

1. Part Number: 1360A-1503

2. Photo:



3. Product Information:

Description: AccuPod® Sample Thief 1500mm, 3 ports

Material of Construction: 316L stainless steel

Method of Construction: Welds ground and polished


Surface Finish: Better than 0.8 microns Ra (tube exteriors)

Overall Length: 1500mm

Outer Tube Diameter: 25mm

Number of Sampling Ports: 3

Position of Sampling Ports:



Nominal Weight: 3.7 Kg

BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin

Recommended Storage Conditions: Dry and ambient temperature

First Use: This sampler should be thoroughly cleaned before being used for the first time

Note: This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24