

**1. Part Number:** 1360A-2001

**2. Photo:**



**3. Product Information:**

**Description:** AccuPod® Sample Thief 2000mm, 1 port  
**Material of Construction:** 316L stainless steel  
**Method of Construction:** Welds ground and polished  
**Surface Finish:** Better than 0.8 microns Ra (tube exteriors)  
**Overall Length:** 2000mm  
**Outer Tube Diameter:** 25mm  
**Number of Sampling Ports:** 1  
**Position of Sampling Ports:**



**Nominal Weight:** 4.9 Kg  
**BSE/TSE statement:** All polishing compounds used in the manufacture of this sampler are of vegetable origin  
**Recommended Storage Conditions:** Dry and ambient temperature  
**First Use:** This sampler should be thoroughly cleaned before being used for the first time  
**Note:** This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

**4. Approval:**

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 14-Nov-24