

**1. Part Number:** 1430S-1500

**2. Photo:**



### 3. Product Information:

Description	Mono Thief 1500mm	
Materials:	Outer Rod:	316L stainless steel
	Tip:	316L stainless steel
	Handle:	316L stainless steel
	Inner Rod:	316L stainless steel
Outer Diameter of Outer Tube:	25 mm	
Overall Length:	1555 mm	
Pocket Length:	135 mm	
Pocket Width:	11 mm	
Pocket Volume:	20 ml	
Number of Pockets	1	
Pocket Position:		



External Finish:	<1 Micron Ra
Internal Pocket Finish:	Fine machined finish
Nominal Weight of Sampler:	5.3 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Nov-24