

**1. Part Number:** 1900A-01

**2. Photo:**



**3. Product Information:**

Description	Measuring Scoop 1ml
Nominal Volume:	1ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Outer Scoop Diameter:	10mm
Inner Scoop Diameter:	7mm
Scoop Depth:	25mm
Overall Length:	155mm
Nominal Weight:	65g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20