


1. Part Number: 1900A-10

2. Photo:



3. Product Information:

| | |
|---------------------------------|---|
| Description | Measuring Scoop 10ml |
| Nominal Volume: | 10ml |
| Material of Construction: | 316L stainless steel |
| Method of Construction: | Crevice free, all welds ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Outer Scoop Diameter: | 25mm |
| Inner Scoop Diameter: | 22mm |
| Scoop Depth: | 30mm |
| Overall Length: | 162mm |
| Nominal Weight: | 90g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature |

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-May-23