

**1. Part Number:** 1900A-10

**2. Photo:**



### 3. Product Information:

Description	Measuring Scoop 10ml
Nominal Volume:	10ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Outer Scoop Diameter:	18mm
Inner Scoop Diameter:	15mm
Scoop Depth:	60mm
Overall Length:	155mm
Nominal Weight:	98g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20