

1. Part Number: 1900A-25**2. Photo:****3. Product Information:**

Description	Measuring Scoop 25ml
Nominal Volume:	25ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Outer Scoop Diameter:	40mm
Inner Scoop Diameter:	34mm
Scoop Depth:	30mm
Overall Length:	185mm
Nominal Weight:	160g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20