


**1. Part Number:** 3044B-100S

**2. Photo:**



### 3. Product Information:

**Description:** Swab Template Kit 100mm

**Kit comprises of:** 1 x Telescopic Rod (1.2 to 3m)  
1 x Elbow  
50 x Swab Templates 100x100mm (sterile)

**Material of Construction:** Telescopic Rods: Aluminium and Plastic  
Elbow: Blue PBT  
Swab Template: Natural HDPE (8888H-100S)

**Template Size:** 100 x 100mm

**Nominal Weight:** 805g

**Recommended Storage Conditions:** Dry and ambient temperature

**Notes:**

- Longer rods are available on request
- All 8888H templates are interchangeable

**Spares:**

Part Number	Description
5065A-3001	Telescopic Rod (1.2 to 3m)
5065G-1B	Angular Elbow (blue)
8888H-50	Swab Template 50x50mm – HDPE (Box of 50)
8888H-50S	Swab Template 50x50mm – HDPE – Sterile (Box of 50)
8888H-100	Swab Template 100x100mm – HDPE (Box of 50)
8888H-100S	Swab Template 100x100mm – HDPE – Sterile (Box of 50)

### 4. Approval:

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 21-Jul-22