

1. Part Number: 3044B-50S

2. Photo:



3. Product Information:

Description: Swab Template Kit 50mm

Kit comprises of: 1 x Telescopic Rod (1.2 to 3m)
1 x Elbow
50 x Swab Templates 50x50mm (sterile)

Material of Construction: Telescopic Rods: Aluminium and Plastic
Elbow: Blue PBT
Swab Template: Natural HDPE (8888H-50S)

Template Size: 50 x 50mm

Nominal Weight: 790g

Recommended Storage Conditions: Dry and ambient temperature

Notes:

- Longer rods are available on request
- All 8888H templates are interchangeable

Spares:

Part Number	Description
5065A-3001	Telescopic Rod (1.2 to 3m)
5065G-1B	Angular Elbow (blue)
8888H-50	Swab Template 50x50mm – HDPE (Box of 50)
8888H-50S	Swab Template 50x50mm – HDPE – Sterile (Box of 50)
8888H-100	Swab Template 100x100mm – HDPE (Box of 50)
8888H-100S	Swab Template 100x100mm – HDPE – Sterile (Box of 50)

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jul-22