

1. Part Number: 3044C-01

2. Photo:



3. Product Information:

Description: Swab Holder Kit

Kit comprises of: 1 x Telescopic Rod (1.2 to 3m)
1 x Elbow
1 x Swab Holder

Material of Construction: Telescopic Rods: Aluminium and Plastic
Elbow: Blue PBT
Swab Holder: 316L stainless steel

Maximum Jaw Opening on Swab Holder: 10mm

Nominal Weight: 990g

Recommended Storage Conditions: Dry and ambient temperature

Notes:

- Swabs are not included with this product
- Longer rods are available on request

Spares:

Part Number	Description
5065A-3001	Telescopic Rod (1.2 to 3m)
5065G-1B	Angular Elbow (blue)
3066A-1	Swab Holder

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jul-22