

**1. Part Number:** 4000A-01

**2. Photograph:**



### 3. Product Information:

Description:	Multilevel Tablet and Capsule Sampler (removeable cups)
Material of Construction:	Frame: 316L stainless Steel Cups: 304 Stainless Steel
Method of Construction:	Cups lift out of frame
Finish	Ra <1.0µm
Number of Sampling Cups:	3
Volume of Sampling Cups:	20ml
Overall Length:	800mm
Nominal Weight:	800g
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 5-Mar-20