

1. Part Number: 4010B-01

2. Photograph:



3. Product Information:

Description:	Multilevel Tablet Sampler
Materials of Construction:	316L stainless steel & PTFE
Method of Construction:	All welds ground and polished
Surface Finish:	Electropolished
Max. Sample Volume:	2.1 Litres
Overall Height:	1275mm
Dimensions of body of Sampler:	1000mm x 95mm x 50mm
Number of Sampling Points:	5
Sampling Heights:	70mm, 270mm, 475mm, 675mm, 880mm
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Weight:	2.8 Kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:


Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 31-Jul-20



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com