

**1. Part Number:** 5005A-100

**2. Photo:**



**3. Product Information:**

|                             |   |
|-----------------------------|---|
| Part No.                    | 5005A-100   |
| Description                 | Stainless Steel Dipper – 100ml  |
| Material of Construction:   | 316L stainless steel  |
| Method of Construction:     | Fully Welded Construction   |
| Surface Finish:             | Cup better than 0.5 microns Ra<br>Handle better than 1.0 microns Ra                     |
| Nominal Cup Volume:         | 100ml   |
| Brimful Cup Volume:         | 100ml   |
| Overall Length:             | 600mm   |
| Max. Width (cup and handle) | 49mm  |
| Cup - Outer Diameter:       | 39mm  |
| Cup - Height:               | 100mm   |
| Nominal Weight:             | 485g  |
| BSE/TSE statement:          | All polishing compounds used in the manufacture of this product are of vegetable origin |

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24