

1. Part Number: 5005A-50

2. Photo:



3. Product Information:

Part No.	5005A-50
Description	Stainless Steel Dipper – 50ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	50ml
Brimful Cup Volume:	50ml
Overall Length:	600mm
Max. Width (cup and handle)	42mm
Cup Diameter:	32mm
Cup Height:	80mm
Nominal Weight:	450g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Feb-20