

1. Part Number: 5005B-500

2. Photo:



3. Product Information:

Part No.	5005B-500
Description	Stainless Steel Dipper – 500ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	500ml
Brimful Cup Volume:	500ml
Overall Length:	1000mm
Max. Width (cup and handle)	85mm
Cup - Outer Diameter:	75mm
Cup - Height:	120mm
Nominal Weight:	980g
BSE/TSE statement: of vegetable origin	All polishing compounds used in the manufacture of this product are

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24