

1. Part Number: 5020L-01S

2. Product Information:

Description: Polyethylene Tubing 1m – Sterile
Tubing Length: 1m (approx.)
Quantity of Tubes per Box: 100
Material of Construction: Low Density Polyethylene (LDPE)
EC Compliance: 1935/2004, 2023/2006, 10/2011, 202/2014
US Compliance: FDA 21 CFR 177.1520
Individually bagged: Yes
Sterile: Yes (Gamma Irradiation 25-45kGy)
Overall Diameter: 8mm
Internal Diameter: 6mm
Pressure Rating: None
ATEX: This product has not been ATEX tested
Box Weight: 3Kg
Recommended Storage Conditions: Dry and ambient temperature

3. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24