

1. Part Number: 5020L-02**2. Product Information:**

Description: Polyethylene Tubing 2m
Tubing Length: 2m (approx.)
Quantity of Tubes per Box 20
Material of Construction: Low Density Polyethylene (LDPE)
EC Compliance 1935/2004, 2023/2006, 10/2011, 202/2014
US Compliance FDA 21 CFR 177.1520
Individually bagged Yes
Sterile N/A
Overall Diameter: 8mm
Internal Diameter: 6mm
Pressure Rating: None
ATEX: This product has not been ATEX tested
Box Weight: 2.2 Kg
Recommended Storage Conditions: Dry and ambient temperature

3. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24