

Product Specification


Specification Number

5020S-01

Revision


02

| | |
|---------------------------------|---|
| Part No. | 5020S-01 |
| Description | Hand Pump Stainless Steel |
| Use: | Reusable liquid pump |
| Material of Construction: | Metallic Parts – 316 stainless steel O-rings – Silicone rubber Handle – Duroplast plastic |
| Sampling Tube Diameter: | OD 8mm |
| Bottle Thread: | The Pump Body is made to suit bottles type 8130A with a 32mm diameter neck* |
| Recommended Storage Conditions: | Dry and ambient temperature |



*Sampling Systems offer a PTFE bottle adapter that can be custom made to suit any bottle with a screw thread

Important Note: If the user continues to operate the pump once the bottle is full the liquid will be drawn into the body of the pump and it will leak out

| | | |
|----------------------------|--|--|
| Specification Approved By: | Signature:  | Name: Jonathan Fenner Position: Director Date: 15-Mar-13 |
| Page 1 of 1 | Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com | |