

**1. Part Number:** 5020T-100USP

**2. Photograph:**



**3. Product Information:**

Description:	PTFE Tubing 100m
Material of Construction:	Polytetrafluoroethylene (PTFE) <ul style="list-style-type: none"><li>• Conforms to: FDA 21 CFR 177.1550</li><li>• Conforms to: USP Class VI</li></ul>
Overall Diameter:	8mm
Internal Diameter:	6mm
Tubing Length:	100m (approx.)
Pressure Rating:	None
Weight:	4.4Kg
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24