

1. Part Numbers: 5030A-1500

2. Photograph:



3. Product Information:

Description	Cup Sampler 100ml - 1500mm
Volume of Cup	100ml (nominal)
Material of Construction:	Cup 316L stainless steel Cup Lid 316L stainless steel Inner Rod 316L stainless steel Outer Rod 316L stainless steel Handles 316L stainless steel Spring Stainless steel (not a product contact part)
Surface Finish:	Satin finish (better than 0.5µm Ra)
Overall Length:	1550mm
Diameter of Cup:	55mm
Weight:	1.3kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24