



	<h1>Sampler Specification</h1>	Specification Number	5035P-600
		Revision	01

Part No. 5035P-600
Description Liquid Master
Material of Construction: Polypropylene (for product contact parts)



Diameter of Outer Tube: 25mm
Max. Diameter of Tip: 33mm
Overall Length: 735mm
BSE/TSE Free: Yes
Nominal Weight: 200g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 23-Apr-13
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	