


**1. Part Number:** 5035T-1001

**2. Photo:**



**3. Product Information:**

Description	Liquid Master PTFE – 1000mm
Materials of Construction:	PTFE & FEP (for product contact parts) <ul style="list-style-type: none"> <li>○ FDA 21 CFR 177.1550</li> <li>○ USP Class VI</li> <li>○ EU regulation 10/2011</li> <li>○ EC regulation 1935/2004</li> </ul>
Outer Diameter of Outer Tube:	25mm
Internal Diameter of Outer Tube:	19mm
Max. Diameter of Tip:	32mm
Diameter of hole in end of tip:	10mm
Overall Length:	1000mm
BSE/TSE Free:	Yes
Nominal Weight:	585g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Oct-19