

1. Part Number: 5050A-24

2. Photograph:



3. Product Information:

Description:	Coliwasa
Material of Construction:	Polyethylene (HDPE) and stainless steel
Method of Construction:	One stainless steel screw is used to attach plug to inner tube One stainless steel screw is used to attach handle to inner tube
Overall Length:	1120mm
Diameter:	21mm
Nominal Weight:	105g
Moulding & Packing Environment:	Assembled & packed in a non-cleanroom factory environment
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	24
Recommended Storage Conditions:	Dry and ambient temperature
ATEX:	This device has not been ATEX tested

IMPORTANT NOTES

- This product is recommended for general industrial use only.
- We do not recommend this product for use in hygienic GMP production areas (food & pharmaceutical).

For food and pharmaceutical use please contact Sampling Systems and ask about:

8050H – SteriWare LiquiThief Sampler

5055H – SteriWare ViscoThief Sampler

8051H – SteriWare Pipettes

SteriWare samplers are made from food grade polymers. They are manufactured and packed under hygienic conditions in a Cleanroom.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jun-22



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com