

1. Part Number: 5065B-1001B

2. Photo:



Angular Beaker fitted to a Telescopic Rod.
Note: Telescopic Rod is sold separately.



Angular Beaker 1000ml
(Blue rotating elbow is included in the price)

3. Product Information:

| | |
|---------------------------------|---|
| Description | Angular Beaker 1000ml – 316Lss |
| Materials: | Beaker: 316L stainless steel Rotating Elbow PBT - Blue |
| Surface Finish: | Better than 0.8 microns Ra |
| Nominal Volume: | 1000ml |
| Internal Beaker Diameter: | 100mm |
| Overall Beaker Diameter: | 110mm |
| Nominal Beaker Height: | 145mm |
| Nominal Weight: | 675g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this sampler are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature |

Note: Suitable for use with 5065A and 5065S telescopic rods.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Dec-24