

1. Part Number: 5065B-1002

2. Photo:



Angular Beaker fitted to a Telescopic Rod.
Note: Telescopic Rod is sold separately.



Angular Beaker 1000ml
(Blue rotating elbow is included in the price)

3. Product Information:

Description	Angular Beaker 1000ml	
Materials of Construction:	Beaker	Polypropylene - Translucent
	Rotating Elbow	PBT – Blue
Nominal Volume:	1000ml	
Brimful Volume:	1100ml	
Maximum Beaker Diameter:	115mm	
Overall Beaker Width:	145mm	
Nominal Beaker Height:	148mm	
Nominal Weight:	181g	
Recommended Storage Conditions:	Dry and ambient temperature	

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Dec-24