

|  |                                |                      |            |
|--|--------------------------------|----------------------|------------|
|  | <h1>Product Specification</h1> | Specification Number | 5065B-2000 |
|  |                                | Revision             | 01         |

Part No. 5065B-2000

Description Angular Beaker

Material of Construction:

|                   |   |
|-------------------|---|
| Beaker            | Polypropylene   |
| Red plastic parts | Polypropylene   |
| Metallic Parts    | Various metals (stainless steel, aluminium, zinc plated mild steel) |




Nominal Volume of Beaker 2000ml

Beaker Diameter: 146mm

Beaker Height: 185mm

Weight: 320g

Recommended Storage Conditions: Dry and ambient temperature

|                            |  |  |
|----------------------------|--|--|
| Specification Approved By: | Signature:<br> | Name: Jonathan Fenner<br>Position: Director<br>Date: 21-Jun-13 |
|----------------------------|--|--|