

**1. Part Number:** 5065B-252

**2. Photo:**



Angular Beaker fitted to a Telescopic Rod.  
Note: Telescopic Rod is sold separately.



Angular Beaker

### 3. Product Information:

Description	Angular Beaker 250ml	
Materials of Construction:	Beaker	Polypropylene - Translucent
	Rotating Elbow	PBT – Blue
Nominal Volume:	250ml	
Brimful Volume:	320ml	
Maximum Beaker Diameter:	73mm	
Overall Beaker Width:	101mm	
Nominal Beaker Height:	109mm	
Nominal Weight:	117g	
Recommended Storage Conditions:	Dry and ambient temperature	

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Jan-23