

1. Part Number: 5065B-252

2. Photo:



Angular Beaker fitted to a Telescopic Rod.
Note: Telescopic Rod is sold separately.



Angular Beaker 250ml
(Blue rotating elbow is included in the price)

3. Product Information:

| | | |
|---------------------------------|-----------------------------|-----------------------------|
| Description | Angular Beaker 250ml | |
| Materials of Construction: | Beaker | Polypropylene - Translucent |
| | Rotating Elbow | PBT – Blue |
| Nominal Volume: | 250ml | |
| Brimful Volume: | 320ml | |
| Maximum Beaker Diameter: | 73mm | |
| Overall Beaker Width: | 101mm | |
| Nominal Beaker Height: | 109mm | |
| Nominal Weight: | 117g | |
| Recommended Storage Conditions: | Dry and ambient temperature | |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 11-Dec-24